Nucleated RBC

Interpretive Summary

Description: Nucleated red blood cells (metarubricytes) represent the early stages of a red blood cell before it matures. They are produced primarily in the bone marrow and are only occasionally observed in peripheral blood.

Decreased Nucleated RBC

Common Causes

• Normal (reference intervals include zero)

Increased Nucleated RBC

Common Causes

- Regenerative anemia
 - Hemolysis
 - Blood loss
- Splenic disease
- Post-splenectomy

Uncommon Causes

- Bone marrow damage
 - Hypoxia
 - Endotoxemia/septicemia
 - $\circ \quad \text{FeLV}$
 - o Drugs
 - o Trauma (fractures)
 - o Myelophthisis: myelodysplastic syndrome, leukemia, metastatic neoplasia
 - Myelofibrosis, bone marrow necrosis
 - Heat stroke
- Extramedullary hematopoeisis, especially splenic
- Lead poisoning
- Iron or copper deficiency
- Hereditary macrocytosis of Poodles

Related Findings

- Hemolysis
 - Increased reticulocytes
 - Increased leukocytes, +/- decreased platelets
 - Increased serum bilirubin, bilirubinuria, +/- hemoglobinuria
 - Spherocytosis (in dogs), autoagglutination, +/- positive Coombs test (IMHA)
 - Positive serology, PCR, or antigen testing for infectious causes
 - Blood parasites visualized on blood smear
 - o Gastrointestinal metallic foreign body found on abdominal radiographs
- Blood loss
 - Increased reticulocytes
 - Decreased total protein and/or albumin
 - o Pleural or peritoneal effusion and/or pulmonary hemorrhage on radiographs or ultrasound
 - Positive fecal ova and parasite screen, positive fecal occult blood



- +/- Decreased serum iron concentration, normal total iron binding capacity, and decreased serum ferritin (if chronic blood loss)
- Increased PT and/or PTT, decreased platelets, prolonged buccal mucosal bleeding time, or low von Willebrand factor level
- Splenic disease
 - o Schistocytes, keratocytes, and acanthocytes
 - Neoplasia or splenitis found on cytology or histopathology
 - o Splenomegaly or splenic mass on abdominal imaging

Additional Information

Physiology

- Nucleated red blood cells are nucleated precursors to mature red blood cells found in circulation.
- Presence in peripheral blood may be associated with a regenerative anemia or due to pathologic conditions affecting bone marrow regulation of red blood cell release.

Diagnostic Methodology

- Circulating nucleated red blood cells are counted along with white blood cells as part of the white blood cell count
 - Automated cell counters cannot differentiate nucleated red blood cells from white blood cells
- Manual slide review is used to count the number of nucleated red blood cells per 100 white blood cells to determine a corrected white blood cell count.
- The correction is performed using this equation: Corrected WBC count = (WBC count [obtained by machine] x 100)/ (# NRBC [obtained during 100 cell differential count] + 100)

References

- Aroch I, Segev G, Loeb E, et al. Peripheral nucleated red blood cells as a prognostic indicator heatstroke in dogs. JVIM 2009; 23(3): 544-51.
- Harvey JW. Atlas of Veterinary Hematology: Blood and Bone Marrow of Domestic Animals. Philadelphia, PA: WB Saunders; 2001.
- Latimer KS, Mahaffey EA, Prasse KW, eds. *Duncan and Prasse's Veterinary Laboratory Medicine: Clinical Pathology*, 4th ed. Ames, IA: Blackwell; 2003.
- Stockham SL, Scott MA. Fundamentals of Veterinary Clinical Pathology, 2nd ed. Ames, IA: Blackwell; 2008.

Last updated 11/1/2013

